



Item	Part Number	Qty	Description
1	DT-082	1	Gallows
2	EG-084/1	1	Head pump
3a	FV-047	1	Metallic handpiece holder Ø 18
3b	EQ-186/A or B	1	Security kit handpiece support
4	EN-049#	1	Surgery handpiece with cord
5	BC-181	1	O'ring Ø 10.00 x 1.00
6	AD-532A	1	Surgery screwed nozzle
7	EK-242	1	1 step foot pedal
8	CD-092	1	Swiss power cable 3 m
	CD-093	1	Hospital grade power cable 3 m
	CD-094	1	European power cable 3 m
	CD-095	1	US-Japan power cable 3 m
9	CE-068	1	Fuses T 1.6 A
10	EQ-015	1	Feet