



MAIA

Sheet	Page
Maia Head	1 - 2
Optical Assembly	3
Central Arm	4
Rear arm with Trasformer	5
Rear arm without transformer	6
Maia Cables	7




Pos	Codice	Description	Validity	pcs
1	SP320901	Optical Cover		1
2	SP320902	Maia Reflector		1 set
3	SP320904	Handle Set		1 set
4	SP320905	Handle Insert Set		1
5	SP320906	Electronic Board	Up to SN:15LD011296	1
5	SP320906A	PCB for lamp with switch	From SN:15LD011297	1
6	SP320926	PCB for lamp with proximity	Up to SN:15LD011296	1
6	SP320926A	PCB for lamp with proximity	From SN:15LD011297	1
7	SP320927	Maia Sensors		1
8	SP320907	Switch Group	Up to SN:15LD011296	1
8	SP320907A	Switch Group	From SN:15LD011297	1
9	SP320909	Cover		1
10	SP320908	Lever Set For Switch		5
11	SP408015	Lable set		3
12	SP320917	Head Cover Set Sensor Type		1
12	SP320911	Head Cover Set		1
13	SP320910	Complete Junction Joke		1


▶ continue on page 2




Pos	Codice	Description	Validity	pcs
14	SP320914	Head Fixing Set		1
15	SP320916	Third Axis Junction Joke		1
16	SP320913	Third Axis Assembly		1




Pos	Codice	Description	Validity	pcs
1	SP320903	Optical Group		1
2	SP320911	Head Cover Set		1
2	SP320917	Head Cover Set Sensor Type		1



Pos	Codice	Description	Validity	pcs
1	SP320951	Central Arm L=550 S		1
1	SP320952	Central Arm L=850 S		1
1	SP320961	Central Arm L=550 TS Unit Version		1
1	SP320962	Central Arm L=550 TS Ceiling Version		1
1	SP320963	Central Arm L=850 TS Unit Version		1
1	SP320964	Central Arm L=850 TS Ceiling Version		1
2	SP320912	Arm Cover Set		1



Pos	Codice	Description	Validity	pcs
1	SP320931	Rear Arm L=750 S Unit Version		1
1	SP320932	Rear Arm L=750 S Ceiling Version		1
1	SP320933	Rear Arm L=900 S Unit Version		1
1	SP320934	Rear Arm L=900 S Ceiling Version		1
2	SP320912	Arm Cover Set		1
3	SP306155	Cover Set for transformer		1
4	SP310405	Transformer		1
5	SP310406	Fuse set		4



Pos	Codice	Description	Validity	pcs
1	SP320941	Rear Arm L=750 TS Unit Version		1
1	SP320942	Rear Arm L=750 TS Ceiling Version		1
1	SP320943	Rear Arm L=900 TS Unit Version		1
	SP320944	Rear Arm L=900 TS Ceiling Version		1
2	SP320912	Arm Cover Set		1



Pos	Codice	Description	Validity	pcs
	SP320075	Cable from Pcb to Remote Control Cable	Up to SN:15LD011296	1
	SP320143	Cable from Pcb to Remote Control Cable	From SN:15LD011297	1
	SP310095	Cable from lamp set to pcb		1
	SP310104	Buzzer		1
	SP320085	Cable from Pcb toPower Supply Cable		1
	SP310139	Cable from pcb to lamp set		1
	SP310146	Remote control cable		1
	SP310075	Power cable for lamp without transformer unit version		1
	SP310403	Power cable for lamp with transformer		1
	SP310413	Power cable for lamp with transformer version		1