


CABLES

page: 8



CENTRAL ARM


see page 5

see page 2

ALYA 2 HEAD


see page 7

REAR ARM
WITHOUT TRANSFORMER




see page 6

REAR ARM
WITH TRANSFORMER







Pos	Code	Description
1	SP310062	Faro lable
2	SP310017006	Led assembly cover white
3	SP310305	Optical group Alya 2nd series
4	SP310025006	Alya 2nd series reflectors
5	SP310140	Insert hook handle set
6	SP310180	Gray handle holder set
6	SP310180006	White handle holder set
6a	SP310352	Grey handle holder set rv2
6a	SP310352006	White handle holder set rv2
7	SP310160	Grey Alya handle set
7	SP310160006	White handle set
7a	SP310370	Grey handle set rv2
7a	SP310370006	White handle set rv2
8	SP310027006	Alya 2nd series complete junction joke
9	SP310029	Grey head joke cover
9	SP310029006	White head joke cover
10	SP310360	Proximity switch with gray cover for Alya second edition
10	SP310360006	Proximity switch with white cover for Alya second edition
12	SP310045	Head fixing set
13	SP310390	Complete handle set rv 2 Gray
13	SP310390006	Complete handle set rv2 White




Pos	Code	Description
1	SP310123	Joystick with gray cover
1	SP310123006	Joystick with white cover
2	SP310102	Joystick Alya
3	SP310108	Lever for joystick switch
4	SP310355	Complete gray Alya switch second version
4	SP310355006	Complete white series switch second version



Pos	Code	Description
2	SP310036006	White cover set
1	SP310451006	Central arm lenght 550mm for Alya-2 with transformer
1	SP310453006	Central arm lenght 855mm for Alya-2 with transformer
1	SP310450006	Central arm lenght 550mm for Alya-2 without transformer
1	SP311450006	Central arm lenght 550mm for Alya-2 without transformer ceiling version
1	SP310452006	Central arm lenght 855mm for Alya-2 without transformer
1	SP311452006	Central arm lenght 855mm for Alya-2 without transformer ceiling version



Pos	Code	Description
2	SP310036006	White cover set
3	SP306155	Edi and Alya transformer cover
4	SP310405	Alya 230V transformer
5	SP310406	Alya 230V fuse set
1	SP310157006	Rear arm lenght 810mm for Alya-2 with transformer
1	SP311401006	Rear arm lenght 810mm for Alya-2 with transformer ceiling version
1	SP310411006	Rear arm lenght 960mm for Alya-2 with transformer
1	SP311411006	Rear arm lenght 960mm for Alya-2 with transformer ceiling version



Pos	Code	Description
2	SP310036006	White cover set
1	SP310150006	Rear arm lenght 810mm for Alya-2 without transformer
1	SP311150006	Rear arm lenght 810mm for ceiling Alya-2 without transformer
1	SP310155006	Rear arm lenght 960mm for Alya-2 without transformer
1	SP311155006	Rear arm lenght 960mm for ceiling Alya-2 without transformer



Pos	Code	Description
1	SP310095	Cable from lamp setto pcb
2	SP310139	Cable from pcb to lamp set
3	SP310138	cable from pcb to power cable
4	SP310104	Alya horn
5	SP310146	Remote control cable
6	SP310075	Alya Ts Power cable
7	SP311075	Alya Ts Power cable ceiling version
8	SP310413	Power cable for lamp S Type